

NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation



NEUROON

OPEN

DOCUMENTATION

v.0.0.1



2017/06/27

in case of technical questions: l.majewski@inteliclinic.com

NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation

Hello!

We are glad you are here! Please enjoy the Neuroon Open.

In case of any technical help or suggestions we will put our best effort to help you.

For more information please visit our community webpage.



Have a nice day,

Inteliclinic team.

IMPORTANT NOTICE

Inteliclinic reserves the right to make corrections, enhancements, modifications improvements and other changes. Documentation may be not error free. Documentation is created for educational purposes. All hardware modifications are strongly not recommended and person who perform them have all responsibility. The aim of this documentation is to improve general knowledge about Neuroon Open hardware what may be necessary for programmers. All product names, logos, and brands are property of their respective owners. All company, product and service names used in this documentation are for identification purposes only. Use of these names, logos, and brands does not imply endorsement.

NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation


Movement Sensor (LIS3DH)

Functionality description:

LIS3DH is MEMS (Micro Electro Mechanical Sensor) ultra-low-power high-performance 3-axis "nano" accelerometer. In NEUROON it is used to detect motion, its duration and intensivity, what is used for wake-up tracking. Remember that head movement usually implies strong muscle signals which may influence the EEG measurement (muscle signals may be even 1000 times stronger than EEG) and sometimes its good to know that disturbances may be a result of movement.

Movement sensor play key role in power management. Using its readings it is possible to put microprocessor to sleep mode and wake it up in most proper time.

For further application it is possible to use it to measure the position of the forehead, however it may be necessary to apply calibration algorithm.



Product note available here:


<http://www.st.com/content/ccc/resource/technical/document/datasheet/3c/ae/50/85/d6/b1/46/fe/CD00274221.pdf/files/CD00274221.pdf/jcr:content/translations/en.CD00274221.pdf>

Features:

- 16-bit data output
 - 12-bit in hi-resolution,
 - 10-bit in Normal resolution,
 - 8-bit in Low-power mode
- I2C interface
- free-fall interrupt generator
- motion detection interrupt generator
- output data rate from 1 Hz to 5.3 kHz
- integrated 32-level FIFO buffer

NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation

Module placement:




LIS3DH module is placed just below Flash memory on the same board as Nordic N51822 SoC.

Hardware connection


| function | pin name | LIS3DH | N51822 |
|----------|-----------|--------|-------------|
| I2C | CS | 8 | VDIG |
| | SDO | 7 | GND |
| | SDA | 6 | 28 |
| | SCL | 4 | 29 |
| | interrupt | INT1 | 17 |

Schematics:



NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation

Temperature Sensors




Neuroon has two independent temperature sensors. Each of them works as independent uncompensated Wheatstone bridge with platinum thermistor PT1000.

PT1000 is placed on bottom layer on separate PCB. PT1000 is very precise and linear thermistor.

Temperature sensors are connected to ADC in NRF SoC

Temperature sensor can be easily simulated in LTSpice like program:




As you may see, in our bridge we used two resistances: 9.1k Ω and 1.1k Ω . It simply means that measurement bridge is in equilibrium when left and right branch has the same resistance. Equilibrium means that the voltage difference between measurement points placed in each branch is equal almost zero. PT1000 has 1.1k Ω resistance in 26 Celsius degrees, what means that our bridge starts to measure temperature when is higher than 26 Celsius degrees.

NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation


Finally the voltage signal (from range 0V up to 3V) enters **10bit Analog to Digital Converter in N51822 SoC**. That simply means that ADC reference has 2 to power 10 ($2^{10}=1024$) voltage levels. Possible references are 1V, 2V, 3V.

| Reference voltage | 1V | 2V | 3V |
|---------------------|------------------------|--------------|-------------|
| Quantisation levels | 1023 (starting from 0) | | |
| Single quant size | $1V/1023=0.0009775V$ | $0.0019550V$ | $0.002932V$ |


Small voltage quants make our thermometer really vulnerable to very small changes in temperature. The highest sensitivity is reached for 1V and 2V reference. Neuroon Open simply changes ADC reference to higher (or lower) depending on reading.




Temperature sensor heat transfer design:



[1]Bare PT1000



[2]Thermal grease



[3]Gold covered electrode

NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation


To avoid allergic skin reactions electrode is covered with thin gold layer. Electrodes are used for both: electrical signal measurements and temperature sensing. To measure the temperature accurately the distance between temperature sensor and electrode is filled with thermal grease.

Functionality description:

Humans are homeotherms, that means they regulate precisely their inner body's temperature. The day-night cycle is regulated by suprachiasmatic nucleus (SCN) located in the same region of the hypothalamus, where "cold" and "hot" reception is integrated, creates a balanced output, like a thermostat. It's no coincidence. Every cell of the mammalian body uses internal circadian oscillators to govern their circadian pattern of gene expression, which is crucial to metabolic activity control. These cellular oscillators are controlled by core body temperature.

Because the temperature might vary due to different reasons, like fever or poor environmental conditions, it is crucial to gather other clinical data before making any assumptions concerning one's health. Changes in core body temperature might be connected to the disturbed circadian rhythm, as observed in high fever, autoimmune diseases. Unfavorable conditions like hot and moist environment might disturb your sleep.



Hardware connections **Amplifier to ADC** displayed below.



NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation

Biological signal amplifier

Biological amplifier circuit plays crucial role. It is responsible for taking measurements of EEG (Electroencefalography) and EMG (Electromyography), due to that fact it need to perform proper signal conditioning. Since the EEG signals are relatively weak (in comparison to muscle signals or noise signals) proper filtration is necessary. Amplification is realised by active, 4-th order low-pass filter. Analogue circuit is placed on separated PCB. Biological amplifier is realised using




Features:

- virtual ground feature
- impedance control circuit (may be used for electrode-skin contact detection and skin hydration sensing)
- EEG adjusted amplifier
- EMG adjusted amplifier
- Notch filter (cancels strong EM components that comes from electrical network)
- active 4-th order lo-pass filter

NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation

- hi pass filter ($f=0.72\text{Hz}$) removes slow electric signals that
- 16bit ADC (ADS1115) 860SPS with internal reference

Schematic diagram




NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation

Pulse oximeter module


Pulse oximeter is realised on separated PCB. It is composed of two LEDs: red and infrared and photodiode. Pulse oximeter sensor is connected to AFE4400 input. AFE4400 is an analog front end for heart rate monitors and low cost pulse oximeters

Features:

- pulse measurement
- monitoring oxygen in blood (usually 92% or less is considered as low level, more than 95% as normal)
- may be used for apnea detection
- integrated LED driver (H-bridge)
- low power
- 13 noise-free bits
- uses SPI Interface



AFE 4400




NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation

Power management and battery charger


Power management is realised using following modules:

- BQ27742 used for battery protection
- two STLQ015M30R used for power distribution
- STC4054GR as battery charger
- STLQ015M30R as LDO
- fuses and thermistors

For user safety batteries are packed into special metal covers. Each battery pack has overheat thermistor attached directly to the cover.




Using USB_CONNECTED pin of NRF you may decide if the device should operate when the USB is attached. We recommend to use USB only for charging and disable normal operation while charging.



NEUROON OPEN : Inteliclinic Fully Nonconfidential Product Documentation

System-on-chip NRF



NRF 51822 is system-on-chip that performs many functionalities. It controls all peripherals, like: battery manager, FLASH memory, LED and vibrator driver, but also realises BlueTooth Low Energy (that's why blue antenna is connected so closely to it). The nRF51822 is built around a 32-bit ARM® Cortex™ M0 CPU with 256kB/128kB flash + 32kB/16kB RAM for improved application performance.

NRF is responsible for saving data on FLASH memory placed just below.

Features:

- 32-bit ARM® Cortex™ M0 CPU with 256kB/128kB flash + 32kB/16kB RAM
- responsible for BLE (Bluetooth Low Energy) communication
- responsible for LEDs and vibrator management
- communicate with external ADC
- its internal ADC is connected with temperature analog front end
- responsible for collecting data from movement sensor

| NRF | function |
|-------------|-----------------|
| SPI | |
| 23 | MOSI |
| 22 | MISO |
| 6 | SCKL |
| I2C | |
| 29 | SCL |
| 28 | SDA |
| UART | |
| 19 | Rx |
| 20 | Tx |

NRF communication interfaces pinout

| NRF | function |
|------------|-----------------|
| AFE | |
| 2 | CS |
| 10 | PDN |
| 30 | Reset |
| 0 | ADC_RDY |

| AFE Diag | |
|-----------------|----------|
| 9 | DIAG_END |
| 8 | LED_ALM |
| 7 | PD_ALM |
| 21 | CLKOUT |

NRF AFE pinout

| NRF | function |
|-----------------------|-----------------|
| Temperature | |
| 4 | AnalogR |
| 5 | AnalogL |
| Electrode test | |
| 3 | Test E1 |
| 1 | Test E2 |
| 14 | Test Gen |

NRF peripheria pinout

| NRF | function |
|---------------------|-----------------|
| Flash | |
| 13 | Flash CS |
| 18 | FlashHold |
| PowerDistrib | |
| 15 | AnalogON |
| 16 | DigitalON |
| other | |
| 25 | LEDS ON |
| 24 | DIRECT ON |
| 11 | USB CONNECTED |
| 17 | ACC INT |
| 12 | POWER BUTTON |